
All-in-One Panel Computer

MLC-S Series

User's Manual

Revision:	1.1
Part No.:	50M-1Z363-1010
Date:	January 8, 2025

This page is intentionally left blank.

1. Contents

1.	Contents	3
1.1	List of Figures	5
2.	Revision History	6
3.	General Information.....	7
3.1	About this Manual	7
3.2	Trademarks.....	8
3.3	Warranty Information	8
3.4	Support, Problems and Failure Analysis.....	8
3.4.1	Recommended action on problems or failures.....	8
3.4.2	Procedure for returning products.....	9
3.4.3	Drivers	9
3.5	Conventions Used in this Manual	9
3.5.1	Conventions for signal words and symbols	9
3.5.2	Categories of persons	10
4.	Safety	11
4.1	Intended Purpose	11
4.1.1	Intended purpose of MLC-S used as stand-alone device	11
4.1.2	Intended purpose of MLC-S combined with compatible devices	12
4.2	Safety and Warning Instructions	13
4.2.1	General electrical hazards	13
4.2.2	Specific hazards caused by electrical leakage currents in devices.....	18
4.2.3	Other hazards	19
4.2.4	Hazards regarding repair and maintenance	26
5.	User Information	29
5.1	General Description	29
5.1.1	MLC-S models	30
5.1.2	MLC-S configurations	30
5.2	Device Overview	32
5.2.1	Interfaces	32
5.2.2	Operational elements	34
5.2.3	Overview front panel buttons.....	35
5.2.4	Functions of the front panel buttons	35
5.2.5	Package contents	39
5.3	Cleaning and Disinfection	40
5.3.1	Safety instructions for cleaning and disinfection of the device.....	41
6.	Specialized Staff Information	43
6.1	Functional Block Diagram.....	43
6.2	Technical Data of Main System Components.....	44
6.3	Detailed System Components	44
6.4	Optional Built-In Components.....	46
6.4.1	Serial ports with galvanic isolation.....	46
6.4.2	WLAN/Bluetooth module.....	46
6.4.3	USB 2.0 ports with galvanic isolation	47
6.4.4	Device configuration with display in portrait orientation (rotation kit)	47

6.4.5	PCIe Riser Card.....	48
6.4.6	RFID	48
6.5	Maintenance and Service	48
6.5.1	General service intervals	48
6.5.2	Initial device installation/cable cover handling	48
6.5.3	Fuse replacement for the internal power supply (AC mains supply).....	51
7.	BIOS Configuration	54
7.1	Default Settings	54
7.2	Enabling PXE/Network Boot	54
7.3	Changing the Boot Sequence	55
8.	Specifications	56
8.1	Mechanical Specifications	56
8.1.1	Front and side dimensions MLC-S	56
8.1.2	Rear dimensions MLC-S	57
8.1.3	Weight.....	57
8.2	Electrical Specifications	58
8.2.1	Mains voltage input.....	58
8.2.2	Power consumption	58
8.3	Environmental Specifications.....	59
8.4	Markings and Labels.....	60
8.4.1	Product identification label.....	60
8.4.2	Mains input fuse label	60
8.4.3	Explanation of label symbols	61
9.	Appendix	62
9.1	EMC requirements	62
9.2	EN 60601-1-2 tables for MLC-S12-21, and MLC-S12-23	65
9.2.1	Guidance and manufacturers declaration – electromagnetic EMISSION	65
9.2.2	Guidance and manufacturers declaration – electromagnetic IMMUNITY	66
9.2.3	Guidance and manufacturers declaration – electromagnetic IMMUNITY – continued	67
9.2.4	Recommended separation distances between portable and mobile RF communications equipment and the MLC-S	68

1.1 List of Figures

Figure 1: MLC-S all-in-one panel computer	30
Figure 2: Bottom view with cable cover mounted	32
Figure 3: Bottom view illustrating device interfaces (cable cover removed)	33
Figure 4: Capacitive touch buttons	34
Figure 5: MLC-S in operation; front touch buttons backlit.....	34
Figure 6: Front touch buttons overview.....	35
Figure 7: Standby status indication.....	36
Figure 8: Power-on status indication.....	36
Figure 9: Display brightness indicator.....	37
Figure 10: Touchscreen disabled indicator	38
Figure 11: All front touch buttons disabled except for the Touchscreen Enable/Disable button	38
Figure 12: Onscreen keyboard activated.....	39
Figure 13: MLC-S series functional block diagram	43
Figure 14: Bottom view illustrating device interfaces	44
Figure 15: Optional components, two COM ports.....	46
Figure 16: Marking label of galvanically isolated USB 2.0 ports with USB symbol and inscription.	47
Figure 17: MLC-S in portrait orientation.....	47
Figure 18: Cable cover installation.....	49
Figure 19: Protective earth strap.....	50
Figure 20: Protective earth strap, I/O interface bezel mounting point.....	50
Figure 21: Protective earth strap, cable cover mounting point.....	50
Figure 22: MLC-S12-xx front and side dimensions in [mm]	56
Figure 23: MLC-S12-xx VESA 100 interface & rear dimensions in [mm].....	57
Figure 24: Product identification label for MLC-S12-XX device	60
Figure 25: Mains input fuse label	60

2. Revision History

Revision	Changes	Date
1.0	Initial release	August 9, 2024
1.1	Address and contact updated	January 8, 2025

3. General Information

3.1 About this Manual

Please read this manual very carefully and assure compliance to all safety instruction and requirements regarding the operation and care of the device. The information contained, will help to prolong the lifetime of the device and to avoid possible hazards.

This manual contains all relevant and necessary information for safe and best possible failure-free operation of the device. However, this manual does not claim to provide gap-free completeness regarding every technical detail or equipment configuration, nor does it include and explain every imaginable installation situation, operating mode or servicing modality. In case any further information not covered by this document is needed, do not hesitate to contact your ADLINK representative.

For further documentation and the latest utilities and drivers, please visit the ADLINK website at www.adlinktech.com or contact your local ADLINK (sales) representative.

The content of this manual shall not become part of or modify any prior or existing agreement, commitment or relationship between ADLINK and the owner/user of the device. The sales contract only describes the entire obligations of ADLINK. The warranty contained in the contract between the parties is the sole warranty of ADLINK.

The information in this document has been checked carefully prior to release. However, no responsibility will be taken for unintended inaccuracies or unintended missing information. Furthermore, ADLINK reserves the right to change any part of this manual without prior notification. ADLINK does not assume any liability for damages, which refer directly or indirectly to the delivery, performance or use of this material.

This documentation is subject to copyright. Reproduction, transmission and use of this document or its contents are not permitted without prior written approval of ADLINK.

Manufacturer

ADLINK Technology Inc

No. 66 Huaya 1st Rd., Guishan Dist.,
Taoyuan City 333411, Taiwan

Email: service@adlinktech.com

Phone: +886-3-216-5088

Internet: www.adlinktech.com

Copyright © 2024, 2025 ADLINK Technology GmbH. All rights reserved. ADLINK will be used as a synonym for ADLINK Technology GmbH in this manual.

3.2 Trademarks

Product names mentioned herein are used for identification purposes only and may be trademarks and/or registered trademarks of their respective companies.

3.3 Warranty Information

Make sure that the device is used exclusively for the intended purpose and is operated and maintained in accordance with the provisions of this manual. Any use beyond the intended purpose is considered non-intended use. The ADLINK product warranty prohibits any modifications to the device.

Your ADLINK product warranty covers material defects and workmanship for the warranty period, starting with the date of shipment. During the warranty period, ADLINK reserves the option to either repair or replace defective products.

Within the warranty period, repair is free of charge if warranty conditions were observed. Due to high expenditures for test procedures, a test fee will be charged if no fault can be found. Repair services outside of the warranty period will be charged.

For further information regarding tests and measurements made in conjunction with warranty and/or repair services, please refer to Section 3.4 below.

3.4 Support, Problems and Failure Analysis

The use of this all-in-one panel computer requires basic IT knowledge, which implies the operation of the system, correct connection of equipment or devices with the available standard interfaces and the handling of the operating system and software.

3.4.1 Recommended action on problems or failures

Each ADLINK product has been tested carefully and thoroughly prior to shipment. In the unlikely event that a problem occurs during operation, please check the BIOS and operating system settings of your system. These settings are often a cause of error due to resource conflicts.

Before contacting ADLINK, please visit the website of ADLINK or the websites of the respective hard- and/or software manufacturers. If the problem cannot be resolved with this documentation and/or updates alone, please contact ADLINK via e-mail or telephone.

Please be prepared to provide the following information:

- Product name, model number and serial number (SN) - see label on the back of the device
- BIOS version - see the version in the BIOS menu or OSD
- Version of the operating system (OS) installed on the device
- Operational environment and system configuration
- Fault description regarding current behavior and expected result
- Report date if the error has been reported already
- The date when the error occurred for the first time

3.4.2 Procedure for returning products

If you need to return your ADLINK device for servicing, please fill in the **Return Material Authorization (RMA)** form before shipping; the RMA form is available on the ADLINK website at https://www.adlinktech.com/en/Medical_Healthcare_Service_and_Support.

Repairs within the warranty period are free of charge if warranty conditions were observed. Due to high expenditures for test procedures, a test fee will be charged if no fault can be found. Repair services outside of the warranty period will be charged.

Your ADLINK product warranty covers material defects and workmanship for the warranty period starting with the date of shipment. During the warranty period, ADLINK reserves the option to either repair or replace defective products.

For warranty service or repair, the product must be returned to a service facility designated by ADLINK, prepaid by the customer.

3.4.3 Drivers

Drivers for your ADLINK device are available at <http://www.adlinktech.com>. Be sure to keep your devices up-to-date with the latest driver and firmware before contacting ADLINK.

3.5 Conventions Used in this Manual

3.5.1 Conventions for signal words and symbols

**DANGER**

Result and cause of the hazard

indicates a hazardous situation which, if not avoided, *will* result in death or serious injury.

**WARNING**

Result and cause of the hazard

indicates a hazardous situation which, if not avoided, *may* result in death or serious injury.

**CAUTION**

Result and cause of the hazard

indicates a hazardous situation which, if not avoided, *could* result in minor or moderate injury and *may* result in damage to property.

**NOTICE**

is used to call special attention to important information related to the device, handling of the device or the respective part of the manual.

3.5.2 Categories of persons

In this manual, the following groups of individuals are categorized:

- Operator
- User
- Specialized staff

Operator

An operator is defined as a natural or legal person, for example, a hospital that owns and has use of the device on whose behalf it is used.

The operator is responsible for providing a safe device and giving users appropriate instructions about the use of the device.

User

Users are professional individuals who meet the training and qualification requirements to operate and work with the device.

It is the user's responsibility to ensure safe use of the device in accordance with the intended use.

Specialized staff

Specialized staff (e.g., an IT administrator employed by the operator or an authorized service technician especially trained by ADLINK) are individuals with adequate training records who have acquired their knowledge through technical education.

4. Safety

ADLINK MLC-S series devices are labelled with the relevant CE mark in accordance with the relevant provisions.

The MLC-S device is manufactured and distributed by ADLINK.

4.1 Intended Purpose

4.1.1 Intended purpose of MLC-S used as stand-alone device

The all-in-one panel computers of the MLC-S series are devices specifically designed to display, monitor and store accumulated data, but not intended to display images on the panel that required accurate color or contrast display, while processing relevant and/or patient data in environments such as:

- surgical OR theatres,
- intensive care units,
- emergency care units, and
- further point-of-care sites

in healthcare facilities such as hospitals, outpatient care units, and other care units.

Handling and operation of the MLC-S series devices shall only be performed by professional personnel of the above-mentioned entities such as medical staff (surgeons, anaesthetists, physicians, etc.), nursing personnel and technical personnel working in healthcare facilities such as IT administrators, engineers, and technical service technicians.

The engineering design and constructional characteristics of the device especially considers requirements for usage within the above-mentioned environment or other hygiene-sensitive environments with respect to:

- hygienic conditions (e.g. fanless, internal cooling concept, screwless housing design),
- quality aspects for displayed / monitored patient data,
- operational control abilities (e.g. with user's hands double gloved),
- processing speed and quality (e.g. endoscopy visualization, etc.), and
- further technical design characteristics (e.g. isolated interface ports, etc.).

The MLC-S series devices are intended for use in indoor environments in accordance with the usage environments described above. The devices are not intended for outdoor use including use in mobile vehicles such as ambulance cars, helicopters, etc. or in other rough environmental conditions of non-domestic environments.

Please refer to the environmental specification defined for safe operation (see Section 8.3).

The described application should not be used for any purpose other than that defined above.

Improper use, alterations made to the devices, or the use of (spare) parts not approved and released by ADLINK may result in unforeseeable damage. Any use of the device not in accordance with the intended purpose will result in loss of all warranty claims. The risks and possible damage resulting thereof are the full responsibility of the user/operator of the device.

4.1.2 Intended purpose of MLC-S combined with compatible devices

Use in conformity with the intended purpose also includes applications of the MLC-S in conjunction with and/or as an integral part of device systems such as anesthesia monitoring systems, respiratory monitoring systems, or various other system applications formed by a combination of the MLC-S with specific software.

Systems which use an MLC-S series device as an integral part or a component are manufactured and placed on the market by third-party companies independently from ADLINK Technology GmbH. The entire system configuration, technical architecture, etc. and the CE conformity procedure for the created system are then the full responsibility of the respective system manufacturer.

Consequently, the manufacturer of the created system determines the explicit intended use of the integrated MLC-S device and/or of the entire system.

4.2 Safety and Warning Instructions

4.2.1 General electrical hazards

**DANGER**

Risk of electrical shock to the user as well as risk of property damage and damage to the device caused by use of inappropriate mains supplies:

Use an original ADLINK mains supply only and preferably the standard internal mains supply only. In case another mains supply must be used, it must follow specifications which provide an equivalent safety level as the original mains supply, as well as technical compatibility to the technical specifications of the device.

Use of a mains supply not authorized by ADLINK may, in addition to the potential hazards caused, result in a loss of warranty!

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**DANGER**

Risk of electrical shock to the user or others and risk of property damage caused by damaged mains connection components:

Never operate the device with a damaged mains supply cable. Damaged mains supply cables including damaged cable plugs must be withdrawn from use immediately. If the mains connection socket of the device is damaged, it must also be withdrawn from operation immediately and must be secured against any further use until the necessary repairs are performed by a trained service technician.

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**WARNING**

Risk of electrical shock or other injuries to the user or others and risk of property damage caused by insufficiently laid mains supply cables:

Take appropriate care in laying out mains supply cables to avoid any possible damage to the cable and to avoid any danger of stumbling over it.

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**WARNING**

Risk of electrical shock or other injuries to the user or others and risk of property damage caused by inappropriate conditions of the mains supply socket:

Never use mains supply sockets which are damaged, dusty, or dirty. Make sure that proper connections are made between the mains connector and the mains outlet socket. Improper plug connections may result in injury caused by sparks and/or a fire hazard.

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**WARNING**

Risk of electrical shock to the user or others and risk of property damage caused by insufficient (ground) protection line connection:

The device must only be operated at mains supply outlets with a protective ground connection. Intentional disruption of the ground protection line and operation of the device with damaged ground protection either in the connecting cable or the device itself are prohibited.

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**WARNING**

Risk of electrical shock to the user or others and risk of property damage caused by overloaded mains supply devices:

Any additional extension cables or multiple outlet sockets being used must be regularly checked to ensure proper condition and secure operation.

Do not overload outlet sockets, extension cables, and multiple outlet socket devices. Doing so may cause electrical shocks or fire hazards.

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**WARNING**

Risk of electrical shock to the user as well as risk of property damage and damage to the device caused by overvoltage:

Appropriate overvoltage protection must be provided to avoid dangerous overvoltage to the device such as those caused by thunderstorms.

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**WARNING**

Risk of electrical shock to the user as well as risk of property damage and damage to the device caused by ingress of fluids or foreign objects:

Although the device is designed with a high level of protection (up to IP54), the device is not fully protected against ingress of fluids. This is the case especially if the device is operated without the protective cable cover at the rear/bottom of the device (over the I/O interface bezel).

To ensure proper protection of the device's housing from fluid ingress, never put any inappropriate objects such as tools, pencils, etc. into openings of the housing (i.e. interface sockets).

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**WARNING**

Risk of electrical shock to the user or others and as well as risk of property damage and damage to the device caused by inaccessible mains sockets:

The device does not provide a mains switch for disconnecting it from the mains supply; the mains disconnection element of the device consists of the mains cable.

Due to this fact, the mains socket and/or the mains cable must always be within easy reach for the user. If the device must be integrated into cabinets or other objects, the operating institution must provide a disconnection switch device at the system level (i.e. operation via an all-lines disconnecting mains switch e.g. at the system).

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**WARNING**

Risk of electrical shock or other injuries to the user or others and risk of property damage caused by insufficient electrical installations:

Electrical installations must be set up by qualified personnel only. For installation and operation of the device, the responsible operating institution must ensure that the chosen personnel is sufficiently qualified with regard to physical, psychological and emotional condition. It is the responsibility of the operating institution to choose qualified personnel for installation and operation of the device.

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4.2.2 Specific hazards caused by electrical leakage currents in devices

WARNING

Risk of an unacceptable high leakage current to the patient, user or others caused by device interfaces not in use:

Never touch unused interfaces such as USB or COM ports and a patient at the same time!

Secure unused interfaces by using appropriate covers, especially the rear/bottom cover delivered with the device and use galvanically isolated interfaces to connect your ADLINK device to other devices.

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WARNING

Risk of an unacceptable high leakage current to the patient caused by insufficient system configuration:

Accessories connected to analog or digital interfaces must be compliant with the respective European or international standards (e.g. EN 62368-1 for data processing devices and EN 60601-1 for medical devices). Furthermore, all configurations for ME systems must be compliant with the EN 60601-1 standard. Users connecting additional devices to a signal input or output are configuring a medical system and are therefore responsible for the system being compliant with the applicable version of the system standard EN 60601-1.

In case of doubt, please contact the technical customer support or your local representative.

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4.2.3 Other hazards

WARNING

Risk of injury to the user or others as well as risk of property damage and damage to the device caused by inappropriate operation conditions:

The device must only be operated within the operation conditions specified by ADLINK.

Never place the device onto heat emitting objects such as radiators or heating ventilators. The environmental temperature must not exceed the maximum temperatures indicated in this documentation. Never place the device into an enclosed housing unless it is ventilated enough to ensure an operational temperature below the allowed maximum.

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WARNING

Risk of injury to the user or others as well as risk of property damage and damage to the device caused by inappropriate operation areas:

The device is **not** intended to be used in surroundings with presence of anaesthetics mixture with air, oxygen or nitrous oxide.

Ensure to install the device locally in an area avoiding these inappropriate surrounding conditions.

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**WARNING**

Risk of injury to the user or others caused by laser emissions in accessory components:

If your ADLINK device is being operated in combination with a product emitting laser light, (e.g. a barcode reader as accessory), never look into the laser beam. This may result in irreversible damage to the eye. Never use settings or functions other than those described in the documentation of the device and never try to repair or disassemble such accessory devices.

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**WARNING**

Risk of injury to the patient caused by inappropriate handling of ceiling mounting equipment:

If the device is mounted to manually adjustable swiveling ceiling supports, always ensure that no collisions with other objects occur when adjustments are made. Harsh collisions may result in damage to the device housing surface by particles that get stripped off thus creating a potential hazard for a patient being present with an open wound.

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**WARNING**

Risk of injury to the user or others caused by substances hazardous to health:

Please note that the device may emit substances hazardous to health (gases, fluids) in case of a fire.

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**WARNING**

Risk of physical injury to the user or others and danger of damage of the device caused by insufficient mechanical mounting:

When mounting the device, use all four VESA 100 threads on the rear of the device with appropriate screws depending on the thickness of the counterpart. Tighten all four screws firmly to support the weight of the device. Make sure that the suspension system can carry the device!

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**CAUTION**

Risk of injury to the user or others, danger of damage to property and danger of damage of the device caused by insufficient mounting or operating places:

The user is responsible for ensuring safe mounting conditions of the device to or on transport devices as well as for observing the safety instructions for the transport vehicles itself. Improper handling may lead to injuries or property damage.

Ensure that the carts, trolleys, surfaces, vehicles, or tables the device is being attached to or onto are stable enough and can bear the applied weight.

Secure the device against accidental movement and from falling.

Please contact your ADLINK partner for further information about appropriate mounting devices.

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**CAUTION**

Risk of injury to the user or others as well as risk of damage to the device caused by lifting and manual transport of the device:

When lifting and transporting the device, consider the product's weight. Move, carry, and put down the device carefully to avoid back or other injuries.

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**CAUTION**

Health hazard to the user caused by allergenic substances:

As with all industrial products, the device may contain substances (e.g. nickel, aluminum) which may induce an allergic reaction. If you suspect an allergic reaction (i.e. skin reaction, frequent sneezing, red eyes, inhalation restrictions or any other similar symptom), immediately consult a physician to determine the cause.

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**CAUTION**

Risk of damage to the device caused by inappropriate usage aids:

If your device is equipped with a touchscreen display, operate this touchscreen only with your finger(s) or an input pen specifically designed for this purpose.

Never use hard, sharp, or tapered objects, e.g. pencils, ball pens, knives, etc. to avoid damage to the touchscreen.

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**CAUTION**

Risk of damage to the device caused by environmental influences:

Never expose your device to excessive vibration, heat or mechanical force. This may lead to damage of your device or parts of it.

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**CAUTION**

Risk of damage to the device caused by aggressive chemical substances:

Although your device is resistant against many aggressive chemical substances, some substances such as acids, chlorine and chlorine dissipating solutions may damage the housing, seals and/or the display.

Please contact your ADLINK representative for further information.

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**CAUTION**

Risk of damage to the device caused by climatic influences:

Never operate your device under condensing conditions which may occur e.g. if the device is brought directly from a cold to a warm surrounding.

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**CAUTION**

Risk of damage to the device itself or devices connected to it caused by influences of electromagnetic fields resulting in the device becoming inoperable:

The device is comprehensively protected against influences of electromagnetic fields by constructional measures (see Section 9.2 regarding electromagnetic immunity). This, however, requires operation of the device only with the cable cover belonging to the device being attached over the I/O interfaces bezel (see Section 5.2.1).

The device is not capable of operating in close proximity to electromagnetic radiation sources which may emit especially high electromagnetic radiation, e.g. certain HF surgery devices, MR devices or cell phones. Ensure a minimum distance of 2 meters from these kinds of devices when operating the MLC-S.

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**CAUTION**

Risk of damage to the device itself or devices connected to it caused by incorrect handling of peripheral devices:

Before connecting your device to additional peripheral devices (e.g. additional monitor screens or other patient monitoring devices), the device must be switched off and disconnected from the mains supply! Disconnect the mains cable from the mains outlet before connecting peripheral devices to your MLC-S computer.

After connection of peripheral devices, the system must only be operated with appropriate potential equalization. Connect the potential equalization pins provided, both the one belonging to the MLC-S, as well as those of the peripheral devices with the potential equalization of the room in which the system is operated in.

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**CAUTION**

Risk of property damage and damage to the device caused by inappropriate system configuration:

Inappropriate changes to the device settings may damage the device or peripheral devices connected to it and may cause the devices or the entire system to function incorrectly. All damage caused by inappropriate settings and/or system configurations will not be covered by warranty!

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**CAUTION**

Risk of property damage respectively environmental damage caused by inappropriate disposal methods:

At the end of the life span of the device it must not be discarded into the general waste, as electronic elements of the device may cause environmental damage if disposed of inappropriately. Please observe the relevant national disposal restrictions valid for your region.

Please contact your local ADLINK partner for further information about device disposal.

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**NOTICE**

Obligation for the user to report incidents:

If the user and/or patient experiences any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**CAUTION**

Risk of damage to the device caused by installation of unapproved function modules:

Installation of function modules via the onboard PCI Express x16, M.2, MiniPCI Express, and 2.5" SATA slots require prior approval by ADLINK and must be installed at an ADLINK production facility. Please contact your ADLINK representative for further details.

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4.2.4 Hazards regarding repair and maintenance**DANGER**

Risk of electrical shock to the user or others as well as risk of property damage and damage to the device caused by dangerous currents:

Hazardous currents are present during operation of the device. Improper operation or maintenance may result in death, severe injuries, or significant damage to property. Strictly observe all instructions being given in this chapter.

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**WARNING**

Risk of electrical shock to the user, patient or others caused by inappropriate repair methods or maintenance:

All repair or maintenance actions must be performed by qualified, ADLINK authorized personnel only! Local respectively national safety regulations and accident prevention regulations must be strictly observed for all servicing activities.

Every change to the device performed by the user or other, not authorized third parties are not allowed and will lead furthermore to loss of device warranty.

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**WARNING**

Risk of electrical shock to the user or others as well as risk of property damage and damage to the device caused by insufficient preparation for repair or maintenance:

Before starting any repair or maintenance procedures, the device must always be disconnected from the mains supply and must be connected to ground.

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**CAUTION**

Risk of electrical shock to the user or others as well as risk of property damage and damage to the device caused by inappropriate spare components:

For any kind of repair or changes to the device (e.g. for post-manufacture upgrades) original spare parts or ADLINK authorized spare parts or components must be used.

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**CAUTION**

Risk of damage to device caused by electrostatic discharge:

The computers, monitors, and their components are highly delicate devices. ESD (electrostatic discharges) and incorrect handling may seriously damage your device.

Only qualified service personnel may handle the device hardware. Connect yourself to ground before opening the device or before touching internal parts!

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5. User Information

All information in this chapter is addressed to the user and operator of the device. It describes all system characteristics, components, and control actions including the maintenance of the device that can be performed by the operator.

Additional information addressed to specialized staff is provided in Chapter 6 “Specialized Staff Information”.

5.1 General Description

The MLC-S series consists of various all-in-one panel computers models which differ in their functional characteristics and, depending on the model, provide different functional elements for operating the devices.

The devices of the MLC-S series are all-in-one panel computers with the enclosure containing the computing elements and the front display combined into a single unit. The operation of computer functions is enabled via the touchscreen interface or via connected peripheral accessories (keyboard, mouse). There are additional functional buttons provided which are described in Section 5.2.2.

The aluminium enclosure is powder-coated with antibacterial paint.

The MLC-S all-in-one panel computer is based on the high-performance 12th generation Intel® Core™ processor family. These devices are best suited for displaying and processing data and operating peripheral devices in applications such as in surgical ORs, intensive care or emergency departments of hospitals and outpatient centers.

The power supply interface of the device is provided by an internal power supply unit,

The VESA 100 interface, located on the back of the device housing, enables mounting of the device onto pedestals or ceiling support units.

5.1.1 MLC-S models

The MLC-S series includes several models. The general product identifier includes a two-digit code, which represents the size of the display diagonal in inch.

Model	Short description
MLC-S12-21	Device with 21.5" Full HD LCD display
MLC-S12-23	Device with 23.8" Full HD LCD display

5.1.2 MLC-S configurations

Figure 1: MLC-S all-in-one panel computer

The MLC-S series is available in several models, all of which have a security glass pane with an optical “anti-glare” coating on the front. This helps minimize unwanted reflections for the user.

The PCAP touchscreen provides easy operation for turning, swiping, zooming and other multi-touch functions like those on modern smart phones (if the software application supports these functions). The touchscreen’s sensitivity is sufficient to provide operation of the touch functionality for a user wearing protective gloves as commonly used in the respective environments.

MLC-S devices are available with different panel sizes (21.5” or 23.8”) and Intel® i-core CPU (i3, i5, i7, i9), see table below. For further information about equipment options that may be provided with your device, refer to Section 6.4 below.

Panel Size		21.5"	23.8"
Computing System	CPU	Intel® Core™ i3-12100TE, 2.10GHz, max 4.00GHz, 4 cores with HTT, 12MB Smart Cache Intel® Core™ i5-12500TE, 1.90GHz, max 4.30GHz, 6 cores with HTT, 18MB Smart Cache Intel® Core™ i7-12700TE, 1.40GHz, max 4.60GHz, 12 cores with HTT, 25MB Smart Cache Intel® Core™ i9-12900TE, 1.10GHz, max 4.80GHz, 16 cores with HTT, 30MB Smart Cache	
	Chipset	Intel® Q670E Chipset	
	Graphics	Intel® UHD Graphics 770 1.45 GHz	
	Memory	Up to 64 GB DDR5 4800 SODIMM in 2 sockets	
Display Characteristics	Resolution	FHD 1920 x 1080	
	Max. Colors	16.7M	
	Viewing Angle	V/H: 89°/89°	
	Touch menu	Optional rotation kit that automatically positions the touch function keys to the bottom of your screen, and your LED bar up when you turn it between landscape and portrait mode.	
	Luminance (cd/m²)	500	<u>625</u>
Storage	SSD	<u>M.2 SSD and / or 2.5" SATA SSD, capacity upon request</u>	
	HDD	<u>Optional, capacity upon request</u>	
I/O Ports	Serial	Up to 6x RS232, 4kV isolated	
	USB	2x USB 2.0 (<u>4 kV isolated optional</u>) 3x USB 3.0	
	DisplayPort	2x DisplayPort, max. resolution 4096 x 2304 @60 Hz	
	WLAN / Bluetooth	<u>Optional WLAN / Bluetooth module</u>	
	LAN	1x 10/100/1000 Mb/s, 4kV isolated 1x 10/100/1000/2500 Mb/s	
	Audio	1x digital microphone	
	RFID	<u>Optional RFID reader</u>	
	Expansion Slot	1x M.2 (2230) slot to support PCIe x1 and USB 2.0 1x M.2 (2280) slot to support PCIe x4 1x miniPCIe (PCIe x1 and USB 2.0) <u>1x PCIe x16 Gen 3.0 slot</u>	
Electricity	Power Supply	100-240V~5A 47-63Hz	

Note: Underlined text indicates selectable or optional configuration.

NOTICE

All of the above listed configurations must be installed at ADLINK's production facilities preferably during the system assembly process. Upgrades after delivery might be possible, please contact ADLINK for further details.

5.2 Device Overview

5.2.1 Interfaces

At the bottom of the device, behind a cover, there are various I/O interfaces allowing connections of the device towards other systems/devices via respective cables and/or connectors. The cable cover protects the interfaces behind it against fluid ingress and, if mounted correctly, provides IP54 ingress protection (against dust ingress and spray water from any angle). The grey plastic elements in the cover allow cables to pass through while maintaining the ingress rating.


Figure 2: Bottom view with cable cover mounted

1. Cable pass-through
2. Cable cover

NOTICE

To operate the device, the cable cover must always be installed to ensure protection against fluid ingress and disturbance by electromagnetic emissions. Please observe the safety instructions in Section 4.2.3.

Underneath the cable cover you will find the device interface ports (see Figure 3 below). These interfaces are identical on all models of the MLC-S series.


Figure 3: Bottom view illustrating device interfaces (cable cover removed)

1. **Fuse bracket** with 2 fuses (5x20mm 2.0A 250V (T), see Section 6.5.3 for details)
2. **Mains connector** (for non-heating device, socket conforms to IEC 60320-1 type C14) for AC mains supply
3. **Speaker**
4. **GbE network interface**, galvanically isolated up to 4 kV (RJ45 connector)
5. **2.5G network interface** (not isolated)
6. **2x USB 2.0** (galvanically isolated version, optional)
7. **Optional: 2x serial ports COM1/COM2**, galvanically isolated up to 4 kV (D-Sub connector 9 pins)
8. **Potential equalization connecting pin** ($\varnothing$ 6mm)
9. **3x USB 3.0** (blue)
10. **2x DisplayPort** (for connection to external displays/monitors)
11. **USB Type-C connector** with Intel® Thunderbolt™ support
12. **Digital microphone**; input interface for external digital microphone
13. **PCIe card expansion slot** (If you'd like to change your configuration, please contact your ADLINK representative for options)
14. **Cover for removable 2.5" SSD** (128 GB standard, other capacities or HDD optional)

NOTICE

(Pos. 8) **Potential equalization connecting pin** ($\varnothing$ 6mm), according to DIN 42801 with color coding yellow/green for connecting a potential equalization cable

The equalization pin of the device must be connected to the potential equalization panel of the respective OR theatre (see IEC 60601-1/EN 60601-1 or the respective national standards).

5.2.2 Operational elements

All MLC-S series devices have a row of capacitive touch buttons at the lower side, under the display's front, with seven function buttons enabling direct control of some device functions.

Figure 4: Capacitive touch buttons

The functional buttons on the front panel of the device are capacitive touch buttons and are activated by the user's finger touching one of the backlit areas. There is no mechanical action with these capacitive touch buttons.


Figure 5: MLC-S in operation; front touch buttons backlit

When the device is powered on, the front panel touch buttons are backlit to indicate the touch-sensitive areas. There are six function buttons backlit in white and the Standby (On/Off) button at the row's center is backlit blue or orange, depending on the operating mode of the device.

5.2.3 Overview front panel buttons


Figure 6: Front touch buttons overview

20. 2x function buttons reserved for future use

21. Call Onscreen Keyboard

22. Standby (On/Off) button for turning on/off the device

23. Touchscreen Enable/Disable button (for cleaning)

24. Display Brightness Decrease button

25. Display Brightness Increase button

When the device is powered on, the front panel touch buttons are backlit to indicate the touch-sensitive areas. There are six function buttons backlit in white and the Standby (On/Off) button at the row's center is backlit blue or orange, depending on the operating mode of the device.

5.2.4 Functions of the front panel buttons

Standby (On/Off) button (22) and status indication

When the device is supplied with power via the mains cable, the Standby (On/Off) button (22) in the center of the front row is illuminated orange to indicate standby status. To turn on the device, touch the Standby (On/Off) button (22) for at least 2 seconds until it illuminates blue and the device will begin the boot routine.

Figure 7: Standby status indication


Figure 8: Power-on status indication

To turn off the device, either:

- use the shutdown procedure of the operating system, or
- touch the Standby (On/Off) button (22) for at least 2 seconds.

In both cases, the device is set into standby mode; the blue indication of the Standby (On/Off) button (22) turns off and changes to orange indication.

NOTICE

While switching off the device via the Standby (On/Off) button (22), do not touch it for more than approx. 7 seconds. Doing so forces a hard shutdown of the device, followed by a restart/boot routine. This may cause a loss of data if it has not been saved or finished processing before shutdown.

NOTICE

To switch off the device completely, it must be disconnected from the mains supply by disconnecting the mains cable. To prevent any data loss, this disconnecting should be done only with the device already in standby mode, i.e. the orange indication is illuminated and the blue indication is off.

Display Brightness Increase (25) and Decrease (24) buttons

If the device is turned on, the brightness of the display can be increased or decreased by touching the Brightness Increase (25) or Decrease (24) buttons continuously or briefly several times until the desired brightness of the display is achieved.

When the display brightness is being adjusted, a green segmented indicator bar is displayed in the lower right corner of the display, on top of the operating system and program application graphics (see figure below). The green indicator bar shows the display brightness setting graphically as well as a percentage value.


Figure 9: Display brightness indicator

To increase the display brightness again, the Brightness Increase button (25) must be activated briefly several times or continuously until the desired display brightness is achieved.

NOTICE

The Display Brightness Decrease button (24) can only decrease the brightness down to 3% of the maximum brightness to avoid the user assuming by mistake the device being switched off when setting it to a very low brightness level. With the display dimmed to the lowest 3% brightness level, minor contour shadows will remain visible.

Touchscreen Enable/Disable button (23)

This function is especially helpful for cleaning and disinfecting the front side of the device while it is in operation. Activating the Touchscreen Enable/Disable button (23) turns the touchscreen off and on again, allowing the device's front panel to be cleaned without interrupting active applications.

To disable the touchscreen function, the Touchscreen Enable/Disable button (23) must be touched for at least 5 seconds until a red symbol appears in the lower right corner of the display (see figure below). This 5-second time delay prevents the user from deactivating the touchscreen unintentionally.


Figure 10: Touchscreen disabled indicator

To reactivate the touchscreen functionality, the Touchscreen Enable/Disable button (23) must be activated for at least 5 seconds until the red symbol in the lower right corner of the display is no longer visible, after which a blue confirmation rectangle appears in the same location for a brief period of time.


Figure 11: All front touch buttons disabled except for the Touchscreen Enable/Disable button

Because of all front panel function buttons are capacitive touch type, all other front panel touch buttons are disabled when the touchscreen is disabled to avoid undesired activation of functions while cleaning the device. The backlight illumination of the other touch buttons (20) to (22) and (24) and (25) are also turned off to indicate that these keys are disabled for the user, until the touchscreen is enabled again.

Call Onscreen Keyboard button (21)

The MLC-S device features an onscreen keyboard for text input via the touchscreen without the need for an external hardware keyboard. The touchscreen keyboard is displayed on top of applications in process.


Figure 12: Onscreen keyboard activated

In the Windows OS, the onscreen keyboard can be called up by swiping the small tab at the left edge of the display onto the screen with a finger. For users wearing a double layer of protective gloves who may find this swipe motion inconvenient or difficult, the Call Onscreen Keyboard button (21) can be activated to display the onscreen keyboard.

Activating the Call Onscreen Keyboard button (21) again will close the keyboard again. Depending on the OS being used, tapping the “close window” function (upper right corner of the window) on the touchscreen of the display is a possible alternative.

Function buttons reserved for future use (20)

These two front panel touch buttons (20) marked as “1” and “2” have no function by default at delivery.

By customer request, these touch buttons can be programmed for specific functions before shipment of the device. Please contact your ADLINK representative for further information (note: possible configurable functions for these keys are highly dependent on the OS and the specific device configuration).

Examples of possible programming of these buttons (if supported by OS and device configuration):

- One button is programmed as a direct “shortcut” to the Windows F1 function, i.e. calls Windows user help.
- One button is programmed to activate a key combination triggering the Windows function for “lock computer”

NOTICE

The customer-specific requirement to individually program the free functional buttons (20) must be specified *before* ordering the device and must be approved by ADLINK’s MLC-S product manager.

5.2.5 Package contents

Please verify that the following item is included in the package:

- MLC-S device

5.3 Cleaning and Disinfection

It is very important that you regularly clean and disinfect your MLC S device with a suitable cleaning or disinfecting agent. The cleaning and disinfection of the device may only be performed by a hygiene specialist or a person instructed by the hygiene specialist. Only ADLINK-tested and approved cleaning agents may be used for cleaning or disinfection. If an agent is not listed in this manual, it should NOT be used as other components could be changed or damaged.

ADLINK recommends periodic product cleaning and disinfection. The following applies to your MLC-S device only. For cables, sensors and other peripheral items, follow the corresponding cleaning instructions of these accessories.

Gently wipe off any dirt on the cabinet or panel surface with a soft cloth soaked in a small amount of deionized water or one of the chemicals listed below. Use a soft, lint-free cloth to wipe the device front until all liquid films are removed. If required, exposure times must be observed for disinfection.

For cleaning the glass / touch front of the device, a simple household glass cleaner available in consumer markets may be used. Use a soft, lint-free cloth to wipe the device front until all liquid films are removed.

Cleaning agents that may be used for disinfection (concentration information per g cleaning solution):

- Glutaraldehyde (max. 2,5 mg/g Glutaral)
- Biguanide (max. 2 mg/g Glucoprotamine)
- Aldehyde (max. 3.75 mg/g Glyoxal / 4.75 mg/g Glutaral)
- Peroxides (max. 4 mg/g Peracetic acid)
- Quaternary ammonium (max. 5 mg/g Benzalkonium chloride / 2.5 mg/g didecyl-dimethylammonium chloride)
- 2-Propanol (60%)

Thoroughly follow the specifications and usage guidelines of the manufacturers of the disinfecting solutions for cleaning your device.

Disinfectants and other agents can lead to material damage of the device, especially if the surface has prior damage (e.g. an accidentally caused scratch) and / or especially over long-term use of these type of disinfectant solutions.

NOTICE

Contaminated items must not be returned without prior written agreement.

5.3.1 Safety instructions for cleaning and disinfection of the device

**WARNING**

Risk of electrical shock to the user caused by presence of potentially conductive liquids:

Prior to any cleaning and disinfection procedure, the mains cable must be disconnected from the mains supply at all times.

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**WARNING**

Risk of unintended operation of the device if equipped with touchscreen:

Before cleaning and/or disinfecting your device, turn off the touchscreen via the Touchscreen Enable/Disable button (23) (see Section 5.2.2 Operational elements).

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**WARNING**

Risk of electrical shock to the user and risk of damage to the device caused by inappropriate processing procedures:

The device is **not** suited for any kind of mechanical and/or thermal decontamination procedures such as autoclaving, steam sterilization, etc.

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**WARNING**

Risk of electrical shock to the user and risk of damage to the device caused by over-dosed amounts of liquids:

Never apply heavy amounts of liquids to your device nor immerse it into liquids of any substance. Be conservative with all cleaning and disinfection solutions used to clean the device and never apply liquids into given openings such as I/O interface ports.

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**CAUTION**

Risk of damage to the device caused by usage of abrasive materials:

Your ADLINK device provides a surface finishing that is sufficiently protected during use within the intended purpose. However, do not use any abrasive materials for cleaning processes, e.g. brushes, sponges, or other objects with hard and/or rough characteristics. Doing so may lead to damage to the front side, the touch surface or the housing surface. Never use e.g. sponges made of metal grafts or similar, nor any sharp or hard objects like knives, pencils etc. for cleaning purposes. Thoroughly follow all instructions provided in this section.

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6. Specialized Staff Information

The information in this chapter is directed to specialized staff (e.g. IT administrator or employees in the IT department) requiring broader technical information beyond the simple operation/handling of the device.

6.1 Functional Block Diagram


Figure 13: MLC-S series functional block diagram

6.2 Technical Data of Main System Components

Component	Technical Data
Processor	Intel® Core™ i3-12100TE, 2.10GHz, max 4.00GHz, 4 cores with HTT, 12MB Smart Cache Intel® Core™ i5-12500TE, 1.90GHz, max 4.30GHz, 6 cores with HTT, 18MB Smart Cache Intel® Core™ i7-12700TE, 1.40GHz, max 4.60GHz, 12 cores with HTT, 25MB Smart Cache Intel® Core™ i9-12900TE, 1.10GHz, max 4.80GHz, 16 cores with HTT, 30MB Smart Cache
Chipset	Intel® Q670E Chipset
Memory	Up to 64 GB DDR5 4800 SODIMM in 2 sockets
Graphics	Intel® UHD Graphics 770 1.45 GHz
Mass Storage	M.2 NVME SSD and/or 2.5" SSD/HDD at SATA 6Gb/s 128 GB SSD is standard equipment (please contact your ADLINK representative for other capacities) 2.5" SSD/HDD drive bay is externally accessible by administrator
Audio Codec	Intel® HD Audio
OS Compatibility	Intel® x64 micro architecture conform

6.3 Detailed System Components


Figure 14: Bottom view illustrating device interfaces

Device Features	
Display	Full HD (1920x1080, 16:9), 16.7M colors, brightness 21.5" with 500 cd/m2 (typ.); 23.8" with 600 cd/m2 (typ.) from panel viewing angle v/h 178°/ 178°
Touchscreen	21.5"/23.8" 10-point, multi-touch Projected Capacitive (PCAP)

Device Features (cont'd)	
Functional buttons, LEDs	Seven PCAP-type functional buttons at the lower side of the front panel consisting of the Standby (On/Off) button with blue/orange backlight and six functional buttons with white backlight (pos. 21-25 in Figure 6) for – display brightness increase/decrease – touchscreen enable/disable (for front cleaning purpose) – calling the onscreen keyboard – two customizable functions
AC mains supply	Internal, wide-range power supply (100...240V~, 50/60Hz); Device protection: fuse bracket with 2 fuses 5x20mm T 2.0A H 250 V, (pos. 1)
Interfaces	
Expansion slot	Slot for PCIe x16 add-on cards (pos. 12) Please contact your ADLINK representative for available options.
Mass storage slot	Covered slot for removable 2.5" SSD/HDD (pos. 14)
Audio	1x digital microphone (pos. 12) 1x speaker (pos. 3)
Display outputs	2x DisplayPort outputs; transferable resolution max. 4096x2304, 60Hz (pos. 10)
USB type-C	USB type-C connector with Intel® Thunderbolt™ support
USB 3.0 ports	3x USB 3.0 SuperSpeed (5Gbit/s) interfaces (blue), compatible with USB 1.1 low (1.5Mbit/s), full (12Mbit/s) and USB 2.0 high speed (480Mbit/s) (pos. 9)
USB 2.0 ports	2x USB 2.0 high speed (480Mbit/s), compatible with USB 1.1 low (1.5Mbit/s) and full speed (12Mbit/s); optional galvanic isolation up to 4 kV (pos. 6)
Network	1x GbE 10/100/1000 Mb/s, 4kV isolated (pos. 4) 1x 2.5G 10/100/1000/2500 Mb/s (pos. 5)
COM serial ports (optional)	COM1 & 2: 2x RS-232 serial ports, max. 256kbps transfer rate; galvanic isolation up to 4 kV (pos. 7) available as optional equipment; on demand, up to 6 serial ports with galvanic isolation are possible as optional equipment if ordered prior to delivery.
Miscellaneous	
Device display configuration	Standard display orientation in landscape orientation; on customer request, display in landscape and portrait orientation available; in both orientation configurations, the functional buttons for directly accessible device functions (see section below) are located at the bottom side of the display and are being automatically activated by rotating the MLC-S device

6.4 Optional Built-In Components

Devices may be ordered with several built-in components options. Retrofit and update of devices already delivered on site is generally not possible. For information about possible retrofit or update, please contact your ADLINK representative.

6.4.1 Serial ports with galvanic isolation

For configurations including two, four, or six COM serial ports, all serial ports have a galvanically isolated design for isolation of up to 4kV.


Figure 15: Optional components, two COM ports

Both serial ports (D-Sub 9 sockets) and the surrounding bracket are mounted into the isolating element also carrying the isolated Ethernet ports.

On demand, up to 3 pairs of 2 RS232 serial ports each with 4kV galvanic isolation are possible as optional equipment, if ordered prior to delivery.

NOTICE

COM 1 and COM 2 cannot be ordered in conjunction with the isolation option of the USB 2.0 ports (see Section 6.4.3). COM 3 and COM 4 cannot be used in conjunction with the PCIe riser card option (see Section 6.4.5).

6.4.2 WLAN/Bluetooth module

If ordered prior to delivery, the devices can be optionally equipped with an IEEE 802.11ac, Dual Band, 2x2 Wi-Fi (PCIe) + Bluetooth 4.0 (USB) module with external antenna(s). The module itself is inserted into an interface on the mainboard.

The used radio frequency bands are 2.4 GHz and 5 GHz depending on settings. Derived from these settings, the bandwidth is 20/40/80 or 160 MHz. Additionally, the transmit power is max. 100mW for the 2.4 GHz band and max. 200mW for the 5 GHz band.

For detailed specifications, please contact your ADLINK representative.

6.4.3 USB 2.0 ports with galvanic isolation

If required, the USB 2.0 ports (pos. 6 in Section 6.3) on the device can be ordered ex-works to provide galvanic isolation up to 4kV. To provide these isolated USB 2.0 ports, an additional circuit board, containing all necessary components is mounted onto the mainboard. These ports cannot be visually distinguished from standard USB 2.0 ports. Therefore, these ports are then marked with a label as shown in the following figure:


Figure 16: Marking label of galvanically isolated USB 2.0 ports with USB symbol and inscription

6.4.4 Device configuration with display in portrait orientation (rotation kit)

If required and ordered in advance, the device can be configured to support the display in portrait orientation in addition to landscape orientation. Then, the front panel function buttons will be automatically activated at the bottom side of the display when rotating the MLC-S.

Figure 17: MLC-S in portrait orientation

6.4.5 PCIe Riser Card

The PCIe x16 slot of the MLC-S CPU board can be used to install a PCIe add-on card in combination with a PCIe riser card. This configuration requires prior approval by ADLINK and must be installed at an ADLINK production facility.

Please contact your ADLINK representative for further details.

6.4.6 RFID

The MLC-S can be equipped with an RFID reader supporting standard ISO 15693 transponders. Please contact your ADLINK representative for further information.

6.5 Maintenance and Service

6.5.1 General service intervals

ADLINK devices do not require servicing or maintenance at regular intervals due to the usage of high-quality industrial standard components. Devices, panel computers, or monitors rarely require servicing when used for their intended purpose.

If your ADLINK system does not operate properly and needs to be returned for servicing, please refer to the instructions in the Sections 3.4.2 and 5.3 above.

In the following Sections 6.5.2 to 6.5.4 you will find instructions for initial device installation and replacing components by specialized staff, if necessary.

6.5.2 Initial device installation/cable cover handling

The device is delivered with the cable cover preinstalled at the rear/bottom of the device. For the initial installation, the mains supply cable and, if required, cables for further external accessory components (e.g. keyboard, mouse), must be fed through the cable cover to connect them to the device.

To connect cables to the device, proceed as follows:

1. Place the MLC-S device with the front side facing down onto a stable working platform protected by a soft layer. Remove the cable cover (Pos.1) at the bottom of the device. The cable cover is fixed by strong magnets located at several points in the bottom side of the device housing. To remove the cable cover, simply grasp it laterally at the housing cutouts left and right, or at the cover's upper edge and pull the cover down from the housing.


Figure 18: Cable cover installation

2. Feed the cables to be connected through the insert pieces (Pos.4) and connect the cables to the respective connectors in the I/O bezel. Please make sure that you choose the suitable insert piece with the appropriate opening diameter for each cable (the bushing set includes a selection of different insert pieces).
3. Insert the pieces with the cables into the cable bushing (Pos.3). Fill a gap which is not used for a cable with an insert without opening.
4. Now push the completed bushing over the mounting plate of the device (Pos.2).
5. Connect the mains supply cable to the mains connector of the device.
6. If the yellow/green protective earth cable inside of the cable cover has been disconnected during installation, reconnect it to the contact provided at the I/O interface bezel of the MLC-S device as shown in Figure 19 to Figure 21 below.
7. Finally, fix the cable cover to the device with “hooking” the cover’s upper edge to the housing first and then tilt the cover towards the housing, until the magnets catch it safely.


Figure 19: Protective earth strap


Figure 20: Protective earth strap, I/O interface bezel mounting point


Figure 21: Protective earth strap, cable cover mounting point

6.5.3 Fuse replacement for the internal power supply (AC mains supply)

**DANGER**

Risk of electrical shock to the user, patient or others as well as risk of property damage and damage to the device caused by hazardous voltage:

Prior to replacing the internal power supply fuses, the device must be disconnected from the mains supply by unplugging the mains cable from the mains outlet.

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**WARNING**

Risk of fire and/or risk of damage to the device caused by inappropriate fuse rating:



The device fuses must only be replaced with fuses of the same rating and type. Carefully observe the device labelling and see also Section 8.4.2.

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
To replace the internal power supply device fuses, follow the steps below:

Step	Illustration	Description
1	N/A	Disconnect the device from the mains supply completely by unplugging the mains connector.
2		Remove the cable cover at the rear/bottom of the device covering the I/O interfaces. Insert a small screw driver or similar object into the notch of the fuse bracket. Pull out the fuse bracket carefully.
3		Turn both fuses towards the sides of the fuse bracket.
4		Pull out the defective fuses.
5		Insert new, identical fuses.

Step	Illustration	Description
6		Push the fuses into place as shown.
7		Insert the fuse bracket and push it into the socket until it snaps into place.
8	N/A	Fix the cable cover to the device.

** NOTICE:**

The fuses of the power supply may blow for different reasons, but this often indicates an internal defect in the device (or a failure of the power supply itself). In case the fuses blow again within a short period after replacement, have your device checked by ADLINK or an authorized ADLINK service technician.

7. BIOS Configuration

7.1 Default Settings

If the BIOS configuration has been set incorrectly and/or the device behaves abnormally, reset the BIOS to its default settings by proceeding as follows:

1. Boot the device and press the key when the POST/boot logo screen appears shortly after the keyboard has been initialized (indicated by flashing keyboard LEDs). The BIOS setup menu will appear on screen.
2. Press the <F9> key
3. Press the <ENTER> key
4. Press the <F10> key
5. Press the <ENTER> key

The device will reboot with the BIOS default settings.

Note: This procedure will also work if the display stays dark and the BIOS itself is not corrupted.

7.2 Enabling PXE / Network Boot

In some cases it may be required to boot the device via a network. To enable the Preboot Execution (PXE) ROM of the integrated Ethernet interface and the correct boot sequence, proceed as follows:

1. Boot the device and press the key when the POST/boot logo screen appears shortly after the keyboard has been initialized (indicated by flashing keyboard LEDs). The BIOS setup menu will appear on screen.
2. Set the cursor to the option **Advanced**
3. Toggle down to the option **Network Stack Configuration** and press the <ENTER> key twice
4. Select **Enabled** and press the <ENTER> key
5. Depending on the local network (IP) requirements, toggle to the appropriate setting and press the <ENTER> key
6. Select **Enabled** and press the <ENTER> key
7. Press the <ESC> key to get back to the main menu
8. Set the cursor to the option **Boot**
9. Set the cursor to the option **Boot Option #1** and press the <ENTER> key
10. Set the cursor to the option **Network** and press the <ENTER> key
11. Set all other listed boot devices to status **Disabled** as described above at point 9
12. Press the <F10> key
13. Press the <ENTER> key to save the settings and to reboot the device

The device will reboot and connect to the network.

Note: Ensure that all necessary services and servers are connected to the network environment.

7.3 Changing the Boot Sequence

As standard setting, the integrated storage medium is the primary boot source. It may be required to boot from another device such as a USB flash drive or external CD/DVD drive to install an operating system (OS) or perform maintenance work. To enable this, the boot sequence can be changed:

1. Connect the required boot source to the device
2. Boot the device and press the key when the POST/boot logo screen appears shortly after the keyboard has been initialized (indicated by flashing keyboard LEDs). The BIOS setup menu will appear on screen.
3. Set the cursor to option **Boot**
4. Set the cursor to option **Boot Option #1** and press the <ENTER> key
5. Choose one of the boot devices listed via the cursor keys and press the <ENTER> key
6. If the boot device is not listed, press the <ENTER> key
7. Set the cursor to option **Hard Drive BBS Priorities** and press the <ENTER> key
8. Set the cursor to option **Boot Option #1** and press the <ENTER> key
9. Choose one of the boot devices listed via the cursor keys and press the <ENTER> key
10. Press the <ESC> key and check if **Boot Option #1** chooses the correct boot source
11. Press the <F10> key
12. Press the <ENTER> key to save the settings and to reboot the device

The device reboots from the previously selected boot source.

Note: If the device cannot boot from the selected boot source, it will attempt to boot from the next boot source in the BIOS start list. To return to booting from the internal storage media, simply disconnect all externally connected boot sources and follow the procedure mentioned above to reinstall the internal storage media as **Boot Option # 1**.

8. Specifications

8.1 Mechanical Specifications

WARNING

Risk of physical injuries and/or property damage:

When mounting the device, use all four VESA 100 threads on the rear of the device with appropriate screws depending on the thickness of the counterpart. Tighten all four screws firmly to support the weight of the device. Make sure that the suspension system is capable of carrying the device!

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8.1.1 Front and side dimensions MLC-S


Figure 22: MLC-S12-xx front and side dimensions in [mm]

Device	a	b	c	d	e
MLC-S12-21	551	342	152	22	85
MLC-S12-23	602	371	166	22	85

Note: All dimensions in [mm].

8.1.2 Rear dimensions MLC-S


Figure 23: MLC-S12-xx VESA 100 interface & rear dimensions in [mm]

Note: For dimension of “c” refer to the previous table.

It is recommended to use M4 screws with a length of 8 mm.

The wall mounting is oriented with the IO side facing down.

8.1.3 Weight

Maximum weight without optional components:

Device version	[kg]	[lbs]
MLC-S12-21	8.9	19.6
MLC-S12-23	10.2	22.5

8.2 Electrical Specifications

8.2.1 Mains voltage input

Input voltage	100-240 V
Input frequency	50-60 Hz
Max. power consumption*	1.5A-0.75A
Fuse (2 pieces):	T 2.0A H 250 V

*Max. power consumption: The documented power values reflect the apparent power measured at the AC input of the PSU.

8.2.2 Power consumption

Test configuration for MLC-S12-23:

Intel® Core™ i7-8700T, clock speed 2.40GHz, 2x 8GB RAM, 128GB M.2 SATA, M.2 2230 frame grabber module, miniPCIe RFID module, 2x galvanically isolated RS-232 ports, 2x galvanically isolated Ethernet ports.

Test configuration for MLC-S-21: as above, but using Intel® Core™ i3-8100T, clock speed 3.10GHz

Model	P_{avg} / P_{max} [W] (load grade 50)*	P_{avg} / P_{max} [W] (load grade 100)**
MLC-S12-21	40 / 44	48 / 50
MLC-S12-23	48 / 72	76 / 80

*load grade 50: Windows® 10 x64 benchmark (10% load for CPU & GPGPU / RAM / 2x SATA / 2D graphics / RFID off), 52% LCD panel brightness; corresponds to a load grade, which represents common average in practical use

**load grade 100: analogous to load grade 50, but respective components under full (100%) load, i.e. Windows® 10 x64 benchmark (100% load for CPU & GPGPU / RAM / 2x SATA / 2D graphics / RFID off), max. LCD panel brightness (100%)

8.3 Environmental Specifications

CAUTION

Risk of device malfunctioning caused by violation of environmental specification limits:

Your ADLINK device is specified for proper operation within the environmental conditions mentioned below. Exceeding these limits can lead to unexpected device behavior and malfunctions.

Always make sure that you operate the device within the environmental conditions specified in this section.

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Ambient temperature range for operation	0°C...+35°C
Ambient temperature range for storage and transport	-10°C...+60°C
Ambient humidity range for operation	10 % ... 90 % (non-condensing)
Ambient humidity range for storage and transport	5 % ... 95 % (non-condensing)
Ambient air pressure for operation, storage and transport	700 hPa ... 1060 hPa

8.4 Markings and Labels

8.4.1 Product identification label

The product identification label is attached to the rear of the device and contains the model and the serial number of the device. The number following “REF” represents the ADLINK product number. The number following “SN” represents the individual serial number of the device.

The following figure shows a generic example of an MLC-S series product identification label. The textual content may differ depending on the respective device characteristics.


Figure 24: Product identification label for MLC-S12-XX device

8.4.2 Mains input fuse label

The label for the mains input fuses is attached at the rear of the device near the mains input socket of the internal power supply.


Figure 25: Mains input fuse label

8.4.3 Explanation of label symbols

Symbol	Meaning
	Read the user's manual and follow all handling and safety instructions.
 	Do not dispose of this device or its parts with conventional domestic waste. Observe your national regulations regarding disposal of electrical and electronic devices. For further information, refer to: DIRECTIVE 2012/19/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 4 July 2012 on waste electrical and electronic equipment (WEEE)
	CE marking: the CE symbol confirms compliance with all relevant guidelines and regulations
 YYYY-MM	Year (YYYY) and month (MM) of manufacture
	Manufacturer
	Catalogue number
	Serial number
IP54	Ingress Protection Code

9. Appendix

9.1 EMC requirements

CAUTION

Adverse events caused by mutual influences between MLC-S and other electrical/electronic equipment:

Although compliant with given EMC standards, the MLC-S might cause adverse events by emission or susceptibility to emissions from other electrical/electronic equipment. Follow the recommendations and information documented in the following EMC-related chapters to avoid or at least mitigate such events. If a mitigation of adverse events is not acceptable, increase the separation distance or relocate the concerned devices.

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NOTICE

To avoid or mitigate unintended emissions and/or irradiations by attached cables use shielded cables with high quality shields. Please observe the recommended bend radius as well. The bend radius is specified by the cable vendor/manufacturer.

CAUTION

Unintended influence of other electrical/electronic equipment due to MLC-S electromagnetic emissions:

The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

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**CAUTION**

Unintended behavior of MLC-S due to emissions of other electrical/electronic equipment such as HF surgical equipment or magnetic resonance imaging equipment:

The operation of equipment or systems with high emission intensity in the vicinity of the MLC-S or attached cables could lead to unintended behavior of the MLC-S. The operation of such devices near the MLC-S is therefore not recommended and is excluded from the electromagnetic environment specified below.

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**CAUTION**

Unintended behavior of MLC-S due to emissions of nearby AM/FM/TV or radar broadcast antennas:

Known RF sources with high emission intensity could cause unexpected or adverse behavior of the MLC-S device. Consider additional shielding measures or contact qualified personnel to solve such issues if relocation or increased separation distance is not sufficient.

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**CAUTION**

Unintended behavior of MLC-S due to emissions of portable RF communications equipment

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the MLC-S, including cables. Otherwise, performance degradation or loss could result.

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**CAUTION**

Device is not tested in a stacked configuration with other equipment:

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

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**NOTICE**

The MLC-S does not need any EMC-related maintenance. All used EMC components such as shields, ground connects and seals are not expected to degrade over the devices lifetime. The tested EM compliance shall be kept therefor.

9.2 EN 60601-1-2 tables for MLC-S12-21, and MLC-S12-23

9.2.1 Guidance and manufacturers declaration – electromagnetic EMISSION

The MLC-S is intended for use in the electromagnetic environment (for professional healthcare) specified below. The customer or the user of the MLC-S should assure that it is used in such an environment.		
Emission test	Device complies to IEC 60601-1-2:2014 EN 60601-1-2:2015 (4th edition)	Electromagnetic environment Guidance (for professional healthcare environment)
RF Emissions CISPR 11	Group 1	The MLC-S uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions Conducted CISPR 11	Class A	The MLC-S is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions Radiated CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliant	

9.2.2 Guidance and manufacturers declaration – electromagnetic IMMUNITY

The MLC-S is intended for use in the electromagnetic environment (for professional healthcare) specified below. The customer or the user of the MLC-S should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2:2014 EN 60601-1-2:2015 (4th ed.) test level	Compliance level	Electromagnetic environment - guidance (for professional healthcare environment)
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transients/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical professional healthcare environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV line(s) to line(s) ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	± 0.5 kV, ± 1 kV line(s) to line(s) ± 0.5 kV, ± 1 kV, ± 2 kV line(s) to earth	Mains power quality should be that of a typical professional healthcare environment.
Power frequency (50, 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	The MLC-S power frequency magnetic fields should be at levels characteristic of a typical location in a typical professional healthcare environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Voltage dips: 0% U_T ; 0.5 cycle 0% U_T ; 1 cycle 70% U_T ; 25/30 cycles Voltage interruptions: 0% U_T ; 250 (50Hz) 300 cycles (60Hz)	Voltage dips: 0% U_T ; 0.5 cycle 0% U_T ; 1 cycle 70% U_T ; 25/30 cycles Voltage interruptions: 0% U_T ; 250 (50Hz) 300 cycles (60Hz)	Mains power quality should be that of a typical professional healthcare environment. If the user of the MLC-S requires continued operation during power mains interruptions, it is recommended that the MLC-S is powered from an uninterruptible power supply or a battery.
NOTE: U_T is the AC mains voltage prior to application of test level.			

9.2.3 Guidance and manufacturers declaration – electromagnetic IMMUNITY – continued

The MLC-S is intended for use in the electromagnetic environment (for professional healthcare) specified below. The customer or the user of the MLC-S should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2:2014 EN 60601-1-2:2015 (4th ed.) test level	Compliance level	Electromagnetic environment - guidance (for professional healthcare environment)
Radiated RF IEC 61000-4-3	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	Recommended separation distance: $d = 1.2 \times \sqrt{P}$ $d = 1.2 \times \sqrt{P} \text{ from 80 MHz to 800 MHz}$ $d = 2.3 \times \sqrt{P} \text{ from 800 MHz to 2.7 GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey [a], should be less than the compliance level in each frequency range [b]. Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted RF IEC 61000-4-6	3 Vrms: 0.15 MHz – 80 MHz outside ISM bands 6 Vrms: in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	3 Vrms: 0.15 MHz – 80 MHz outside ISM bands 6 Vrms: in ISM bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the MLC-S including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
[a] Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MLC-S is used exceeds the applicable RF compliance level above, the MLC-S should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the MLC-S. [b] Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

9.2.4 Recommended separation distances between portable and mobile RF communications equipment and the MLC-S

The MLC-S is intended for use in an electromagnetic environment (for professional healthcare) in which radiated RF disturbances are controlled. The user of the MLC-S can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MLC-S as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = 1.2 \times \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \times \sqrt{P}$	800 MHz to 2.7GHz $d = 2.3 \times \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.79	3.79	7.27
100	12.00	12.00	23.00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters [m] can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.