




Open eVision Deep Learning Studio

Deep Learning training and evaluation application

At a Glance



- Ease the evaluation of Open eVision's Deep Learning tools
- Dataset creation and image annotation for classification, segmentation, and object localization
- Create and configure dataset splits to decide how your images are used
- Manage the data augmentation transformations
- Train your tools in succession thanks to the training queue
- Validation and analysis of the results of the trained tools
- Available on Windows and Linux
- Free of charge

Benefits

New in Open eVision 23.08

- Barcode grading with EasyBarCode2
- Benchmarking tool in Deep Learning Studio
- Learn feature for EasyMatrixCode
- Detect damaged data matrix codes
- Grading and better decoding for EasyBarCode
- Speed improvements for EasyMatch

Open eVision Deep Learning Studio

Open eVision Deep Learning Studio is an application that assists the user during the creation of the dataset as well as the training and testing of the Deep Learning tools of Open eVision.

Open eVision Deep Learning Studio is free and does not require any license. It allows you to test the Deep Learning libraries using your own images. No programming is required.

Just click on **DOWNLOAD OPEN EVISION DEEP LEARNING STUDIO** and install Open eVision. Sample images, manuals and example programs are included.

Deep Learning projects

A Deep Learning Studio project manages your dataset and the Deep Learning tools you created. A project is associated with one the Deep Learning tool (EasyClassify, EasySegment Unsupervised, EasySegment Supervised or EasyLocate) and supports all their features.

Within a project, you can create as many tools as you want. It allows you to easily experiment with the parameters of the tools, different split of the dataset, or data augmentation settings.

New in Open eVision 23.04

- A new code reader tool ECodeReader allowing to read multiple code types from a single integrated interface
- Deep Learning Studio is now available on Linux ARM-64 platforms
- Neural Network “engines” allow additional devices support (OpenVINO, Tensor RT) and increased speed.
- EasyBarCode2 now supports learning from images to improve its capabilities on similar images.

Control how your images are used

Deep Learning Studio allows you to split your dataset into training, validation, and test sets. You can create multiple dataset splits to experiment and check the performance of tools trained with different set of images.

You can create dataset splits at random or manually select the set of each image.

Data augmentation

The rich data augmentation capabilities of EasyClassify, EasySegment, and EasyLocate are available in Deep Learning Studio. Tune and visualize the geometric, color, and noise data augmentations. You can create different set of data augmentation settings to experiment how it influence your results.

Configure and train your tools

The Tools tab allows you to configure and train your tools. Operating on CPU or GPU, the training can be stopped and restarted at any time. You can launch as many training as you want thanks to the processing queue. The training and inference operations will be queued and processed one after the other.

Validation and result analysis

The validation process is customized for each library to allow you to get the most out of your data. A comprehensive set of metrics, tables, and/or graphs is available to analyze and explore the results of the training process.

Tables and confusion matrixes allow you to filter your results to understand the strengths and weaknesses of the trained models. Score histograms and ROI curves are useful to select a threshold and adapt the trained models to your needs.

Annotate your dataset

Deep Learning Studio integrates annotation tools adapted to each library. For classification and unsupervised segmentation, you can quickly assign label to each image. For supervised segmentation, the segmentation editor allows you to draw the ground truth segmentation. For localisation, the object editor allows you to quickly draw the bounding box around each of your objects.

The image editor also allows you to select a region of interest and mask parts of your image.

New Assisted Segmentation Tool

- Assisted segmentation tool to ease and speed up annotation.

Specifications

Software

Host PC Operating System	• Microsoft Windows 11, 10, 8.1, 7 for x86-64 (64-bit) processor architecture • Microsoft Windows 10, 8.1, 7 for x86 (32-bit) processor architecture • Linux for x86-64 (64-bit) processor architectures • Minimum requirements: – 8 GB RAM – 400 MB free hard disk space
--------------------------	--

Ordering Information

Product code - Description	• 4912 - Open eVision Deep Learning Studio
----------------------------	--



EMEA

Euresys SA

Liège Science Park - Rue du Bois Saint-Jean, 20
4102 Seraing - Belgium

Email: sales.europe@euresys.com

EMEA

Sensor to Image GmbH

Lechtorstrasse 20
86956 Schongau - Germany

Email: sales.europe@euresys.com

AMERICA

Euresys Inc.

316 Prado Way
Greenville, SC 29607 - United States

Email: sales.americas@euresys.com

ASIA

Euresys Pte. Ltd.

750A Chai Chee Road - #07-15 ESR BizPark @ Chai Chee
Singapore 469001 - Singapore

Email: sales.asia@euresys.com

CHINA

Euresys Shanghai Liaison Office

Unit 802, Tower B, Greenland The Center - No.500 Yunjin Road, Xuhui District
200232 Shanghai - China

Euresys上海联络处

上海市徐汇区云锦路500号绿地汇中心B座802室
200232

Email: sales.china@euresys.com

CHINA

Euresys Shenzhen Liaison Office

Room 1202 - Chinese Overseas Scholars Venture Building
518057 Shenzhen - China

Euresys深圳联络处

深圳南山区留学生创业大厦1期1202
518057

Email: sales.china@euresys.com

JAPAN

Euresys Japan K.K.

Expert Office Shinyokohama - Nisso Dai 18 Building, Shinyokohama 3-7-18, Kohoku
Yokohama 222-0033 - Japan
〒222-0033

神奈川県横浜市港北区新横浜3-7-18 日総第18ビル エキスパートオフィス新横浜

Email: sales.japan@euresys.com

More at www.euresys.com

